

Sabina™ II EE sit-to-stand lift

Versatile design with enhanced adaptability



The Sabina II EE sit-to-stand lift is designed for mobilising patients not yet ready to fully stand on their own.

Features and benefits

Caters for a variety of patients with **two different slingbar options available**, as well as a selection of sit-to-stand vests

Optional comfort slingbar, combined with the Liko ComfortVest provides an especially gentle lifting action

An extensive assortment of accessories is available including side support, slingbars, heel support and seatstrap

Adjustable height settings accommodate a variety of different patient sizes

Powerful lifting motor ensures the smooth raising and lowering of patients

Base-width adjustment offers greater flexibility of use when lifting from a bed or chair, enhancing patient mobility

Sabina II EE sit-to-stand lift



Height adjustments



Heel cups for enhanced patient safety



Assortment of Liko vests

Technical Specifications

Maximum Load

active lifting	200 kg
passive lifting	150 kg

Weight	45.2 kg
---------------	---------

Protection class	IPX4
-------------------------	------

Slingbar dimensions	350mm Slingbar or Comfort 600mm slingbar
----------------------------	---

Base

height	107 mm
length	1060 mm
width	690 - 1115 mm
expansion	electric

Castors

front	75 mm double band
rear	75 mm single band

Material	powder coated steel
-----------------	---------------------

Emergency lifting	electric
--------------------------	----------

Emergency lowering	manual and electric
---------------------------	---------------------

Lifting speed (without load)

with Slingbar 350	54 mm / sec
with Comfort slingbar	47 mm/sec

Lifting interval

with Slingbar 350	825 mm
Comfort slingbar	785 mm

Sound level	46 dBA
--------------------	--------

Battery	lead-acid gel
----------------	---------------

Charger	built-in
----------------	----------

Hand control	wired
---------------------	-------

Product number	2020003
-----------------------	---------

Hill-Rom is a leading global medical technology company with approximately 10,000 employees worldwide. We partner with health care providers in more than 100 countries by focusing on patient care solutions that improve clinical and economic outcomes in five core areas: **Advancing Mobility, Wound Care and Prevention, Patient Monitoring and Diagnostics, Surgical Safety and Efficiency, and Respiratory Health.** Around the world the people, programs, and products of Hill-Rom work towards one mission: **Enhancing outcomes for patients and their caregivers.**

Not all products/options are available in all countries. For further information about our products and services, please contact your local Hill-Rom representative or visit our webpage:

www.hill-rom.com

Sabina II EE mobile lift is intended to be used for the lifting and transfer of patients in a variety of care settings.

Class I

Manufacturer: Liko AB, Nedre vägen 100, SE-975 92 Luleå, Sweden

This medical device is a regulated health product which, pursuant to such regulation bears a CE mark. Hill-Rom recommends that you carefully read the detailed instructions for safe and proper use included in the documents accompanying the medical devices. The personnel of healthcare establishments are responsible for the proper use and maintenance of these medical devices.

Hill-Rom reserves the right to make changes without notice in design, specifications and models. The only warranty Hill-Rom makes is the express written warranty extended on the sale or rental of its products.

©2017 Hill-Rom Services, Inc. ALL RIGHTS RESERVED.

Doc. No: 5EN183314-01, 09 February 2017

Enhancing outcomes for patients and their caregivers:

Hill-Rom